

Thermo-transfer printer indicator labels

Documentation



It sticks .. and conform to standards

Almost every Thermo-transfer printer indicator label is a customised product, manufactured to each customer's specific needs.

Options for Thermo-transfer printer indicator labels:

- ★ Single or double (piggy back) label material
- ★ Tyvek, PP or Vellum top material
- ★ Indicators for STEAM, EO, FORM, Plasma
- ★ Colored labels (incl. logo or Hospital address data)
- ★ One, two or three part label
- ★ On roll, folded or as sheets

The Thermo-transfer printer labels are especially selected from raw material suppliers to suit the harsh environment conditions of steam sterilization. The selection was based upon requirements for the glue, the top material and ink printing conditions.

After serious and severe testing we found the right raw materials which have been tested by us both under laboratory conditions as well as out in the field.




Majac Medical Products Pty Ltd
Ph: +617 3265 6355
Fax: +617 3865 2729
sales@majacmedical.com.au
www.majacmedical.com.au



FS 540017

Also Available From:



Thermo-transfer printer indicator labels

Documentation



Any requirement we can fulfill

The Steritech™ Label material testing was based upon the following requirements:

- ★ Stickiness, even on difficult adhering surfaces
- ★ Behavior under high temperature and relative humidity
- ★ Behavior of printed ink under these conditions
- ★ Natural Rubber Latex free
- ★ REACH requirements
- ★ Free of any toxic particles
- ★ Suitable for printing by Thermo-transfer printers
- ★ Label gun printing ink acceptance and drying time

The Steritech Label materials fulfill and meet most of our customers requirements and now and then supersedes these.



Majac Medical Products Pty Ltd
Ph: +617 3265 6355
Fax: +617 3865 2729
sales@majacmedical.com.au
www.majacmedical.com.au



Also Available From:

