

CABINET COUCH WITH SLIDE-OUT STEP STOOL: 2 SECTION

Code: 1205

- Cabinet made from colour-board with impact resistant PVC edging
- Sturdy construction and quality finish
- Cabinet colours: lemon tree, beech, jarrah, grey
- Vinyl upholstery colours: Windsor blue, charcoal, black or doeskin
- Custom cabinet and vinyl colours available
- Five deep drawers
- Integrated pull-out step stool
- Retractable instrument shelf
- Adjustable headrest: 580mmL x 585mmW
- Overall dimensions: 1790mmL x 580mmW x 830mmH
- Safe working load 150kg
- Matching single and double step-up stools available

