


Cabinet Couches

- 1830mmL x 585mmW x 820mmH
- Safe working load 150kg
- Cabinet made from colourboard with impact resistant pvc edging
- Sturdy construction and quality finish
- Sliding panel doors with inner shelf
- Five drawers
- Retractable instrument shelf
- Custom colours available



Model 1203 Combination Timber Couch

Cabinet:



Vinyl:



Other colour combinations available

Cabinet Step-Up Stools

Model 1220 Single Step-Up Stool

Model 1221 Double Step-Up Stool

Timber colour to match cabinet couch

