

“GP” Examination Bed - Height Adjustable

2 Section Exam Couch

“GP” Examination Couch (Table) – Height Adjustable, Designed for General Practice

Meets the minimum specifications of disability advocates (www.racgp.org.au/standards/513)

TGA approved ARTG No. 123573

Code: 20.07121 with manual backrest

Code: 20.07121.G2 with gas strut backrest



The “GP” Exam Couch - (electric height and manual backrest). This model is a versatile examination couch.

Specifications:

- * Grey base / Wedgewood upholstery / 710mm wide
- * Overall length 1930mm / backrest section 680mm [0/+80°]
- * Electric height adjustment [45 - 95cm] with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200 kg

Grey Base / Wedgewood upholstery. Condor brand medical grade heavy-duty vinyl fabric, delivers excellent strength and wear characteristics.

Formulated to inhibit hardening and breakdown of the surface due to perspiration and oils; fire retardant; sanitised; stain resistant.



Manual back rest. Breathing Hole and plug is a factory fitted option for additional \$75.00. TGA approved ARTG No.123573.

Note: Breathing Hole optional

Hospital Examination Couch

Hospital Examination Couch – 2 Section

Code: 20.07121H with manual backrest



- The Hospital Examination Couch has been developed in Australia to satisfy the requirements of the RACGP recommendations.
- As such, it is long enough and wide enough for most patients. The minimum height of just 45cm makes access for frail or disabled patients easy and the maximum height of 95cm is high enough for tallest carers to work comfortably. Accessory cot sides and gynaecological stirrups can also be fitted.
- The 125mm braked castors allow easy movement and 155mm clearance below the frame for patient hoist access.
- The height and backrest are electrically adjusted using a hand control.

Options:

Upholstery Width – 800 mm wide +\$120.00
Breathe Holes – in Breathe Hole in backrest with plug +\$75.00

Upholstery – Leather +\$400.00

Mobility Options – High Retractable castors (155mm clearance for hoists) no extra charge

Lifting capacity – 350kg Lifting +\$605.00

Upholstery angle adjustment

- Electric Backrest (with hand control) +\$205.00
- Electric Backrest (with foot control) +\$330.00
- Gas strut assistance on Backrest +\$205.00

Electric Control – Twin foot or hand control +\$325.00

Battery – Add Battery +\$455.00

Upholstery Colour Options:

(Colours are approximate only)

Studio Encore Standard Vinyl Colours



Snow




Chalk



Brazil



Lemon Twist



Lime



Tangerine



Flame



Seagrass




Caribbean



Rhubarb




Lipstick




Iceburg



Wedgewood



Charcoal



Smoke



Black

Colours are approximate only depending on your monitor

Condor brand medical grade heavy-duty vinyl fabric. The exclusive formulation delivers a superior softness for comfort including excellent strength and wear characteristics. Formulated to inhibit hardening and breakdown of the PVC surface due to perspiration and oils; light fast; fire retardant; sanitised; stain resistant.

Maintenance

Recommended cleaning with chemically neutral soap and warm water, then wipe clean. Do not dry-clean or use abrasive or solvent based products.