

“GP” Examination Bed - Height Adjustable

3 Section Exam Couch

“GP” Height Adjustable Examination Couch / Treatment Chair – 3 Section

Meets the minimum specifications of disability advocates (www.racgp.org.au/standards/513)

TGA approved ARTG No.123573 Code:30.07121



Examination Couch – Height Adjustable – 3 Section.

Specifications:

- * Grey base / Wedgewood upholstery / 710mm wide
- * Three sections: 360mm with breathe hole & plug (+50 / -80°) / 560mm / 940mm (+80deg) – Overall length 1930 mm
- * Electric height adjustment (45 – 95 cm) with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200kg
- * Breathing Hole and plug is a factory fitted option for additional \$75.00.

Note: Breathing Hole optional

This versatile chair can also be used as an examination bed


Height adjustable Arm Rests

Optional on Manipulation and examination couches.
Converts a three section couch for multi-purpose use.
Can be retrofitted on-site to suitable couches

As the backrest lifts, the seat tilts backwards to provide a more comfortable configuration.

Hospital Examination Couch – 3 Section

Code: 30.07121H



- The 3 section Hospital Examination Couch has been developed in Australia to satisfy the requirements of the RACGP recommendations.
- As such, it is long enough and wide enough for most patients. The minimum height of just 45cm makes access for frail or disabled patients easy and the maximum height of 95cm is high enough for tallest carers to work comfortably. Accessory cot sides and gynaecological stirrups can also be fitted.
- The three ergonomically sized sections allow for more flexibility in use from a simple plinth through to a chair configuration.
- The 125mm braked castors allow easy movement and 155mm clearance below the frame for patient hoist access.
- The height and head rest and long section are electrically adjusted using a hand control.

Options:

Upholstery Width – 800 mm wide +\$120.00

Breathe Holes

- Breathe Hole with plug in Short Section +\$75.00
- Breathe Hole with plug in Long Section +\$75.00
- Breathe Hole with plus in both Long and Short Sections +\$150.00

Upholstery – Leather +\$400.00

Mobility Options – High Retractable castors (155mm clearance for hoists) no extra charge

Lifting capacity – 350kg Lifting +\$605.00

Electric Control – Foot control +\$415.00

Battery – Add Battery +\$455.00

Upholstery Colour Options:

(Colours are approximate only)

Studio Encore Standard Vinyl Colours



Snow




Chalk



Brazil



Lemon Twist



Lime



Tangerine



Flame



Seagrass




Caribbean



Rhubarb




Lipstick




Iceburg



Wedgewood



Charcoal



Smoke



Black

Colours are approximate only depending on your monitor

Condor brand medical grade heavy-duty vinyl fabric. The exclusive formulation delivers a superior softness for comfort including excellent strength and wear characteristics. Formulated to inhibit hardening and breakdown of the PVC surface due to perspiration and oils; light fast; fire retardant; sanitised; stain resistant.

Maintenance

Recommended cleaning with chemically neutral soap and warm water, then wipe clean. Do not dry-clean or use abrasive or solvent based products.