EXAM COUCH: 2 SECTION

EXAM COUCH WITH POWDERCOATED FRAME

CODE: 1000

- STRONG POWDERCOATED STEEL FRAME
- VINYL UPHOLSTERY COLOURS: GREY OR BLACK
- FLAT PACKED
- SIMPLE AND EASY ASSEMBLY, 4 BOLTS AND HEX TOOL PROVIDED
- ADJUSTABLE HEADREST: 610MML x 610MMW
- OVERALL DIMENSIONS: 1850mmL x 610mmW x 800mmH
- SAFE WORKING LOAD 175KG





Majac Medical Products Pty Ltd

Ph: +617 3265 6355 Fax: +617 3865 2729 sales@majacmedical.com.au www.majacmedical.com.au

