

Examination Couches & Step Up-Stools



Standard Exam Couch

- Safe working load 125kg
- Adjustable headrest
- Flat packed for self assembly
- 1820mmL x 610mmW x 780mmH

Model 1001 Standard Exam Couch Powdercoated

Vinyl:     Frame:  

Model 1002 Standard Exam Couch Chrome Plated

Vinyl top:  



Heavy Duty Exam Couch - 250kg SWL

- Safe working load 250kg
- Comfortable thick foam padding
- Adjustable headrest
- Flat packed. Tool-less assembly
- 1830mmL x 610mmW x 800mmH

Model 1004 Heavy Duty Exam Couch

Vinyl:   On grey frame only

Model 1048 Drawers set of 3 for 1004

Optional bolt-on metal sided drawer set



Step-Up Stools

- Colour combinations to match couches
- Non slip vinyl tops with rubber edges
- Single size 455mmL x 265mmD x 200mmH
- Double size 455mmL x 550mmD x 470mmH

Model 1060 Single Step-Up Stool – Chrome

Model 1061 Single Step-Up Stool – Powdercoated

Model 1061HR Single Step-Up Stool – With Hand Rail

Model 1062 Double Step-Up Stool – Chrome

Model 1063 Double Step-Up Stool – Powdercoated

Model 1063HR Double Step-Up Stool – With Hand Rail