

**LIFE+GUARD™**

**WEARABLE TOUCH-FREE  
SANITIZER DISPENSER**



Sanitize Your Hands Anytime + Anywhere

**FOR MEDICAL PROFESSIONALS**



## Innovative wearable hand sanitizer dispenser

- ✓ Hand sanitizer always within reach “ANYTIME+ANYWHERE”
- ✓ Worn on your belt or on a lanyard
- ✓ Touch-free operation. Simply place 2 fingers at the bottom and an infra-red beam activates a built-in pump that dispenses sanitizer into your hand.
- ✓ Refill with hospital preferred hand sanitizer gel
- ✓ Contains enough sanitizer for a whole day (60ml capacity)
- ✓ Single recharge delivers 1,000 doses
- ✓ Easy recharge with a USB cable (supplied)
- ✓ 3 electronic dispense volume settings
- ✓ Sanitizer level window
- ✓ Electronic battery status indicator



## A WEARABLE DISPENSER PUTS SANITIZER ALWAYS WITHIN REACH

LIFE+GUARD is designed to be worn on the body and has belt and lanyard mounting brackets included.

LIFE+GUARD is never out of reach and is accessible at all times to the healthcare professional.

Easy access to sanitizer promotes good hand hygiene, which benefits both patients and the staff.

This wearable device was designed for medical professionals on the move. Easy to refill, easy to recharge.




## FIGHT AGAINST HOSPITAL ACQUIRED INFECTIONS

Practicing hand hygiene is a simple yet effective way to prevent infections.

On average, healthcare providers clean their hands less than half of the times they should.

On any given day, about one in 31 hospital patients has at least one healthcare-associated infection.

LIFE+GUARD puts hand sanitizer constantly within reach, making it simple for healthcare workers to clean their hands when they need to.



World Health Organization's five moments for hand hygiene

### **Professional Applications**



Hospital applications include:

- Medical staff on ward rounds, moving quickly from patient to patient
- Nurses giving patient care
- ‘Hands-on’ applications such as physiotherapy
- Elderly care facilities
- Hospital Porters
- Hospital cleaning and maintenance staff
- Food service workers serving patient meals

**NO INSTALLATION OR MAINTENANCE  
REQUIRED, UNLIKE FIXED DISPENSING  
SYSTEMS**

## LIFE+GUARD PROFESSIONAL USERS

### MOBILE

Always on the go



### SHARED SURFACES

In contact with objects touched by others



### PHYSICAL CONTACT

Interacting with patients and staff



Keep your staff and patients safe.

Utilize your approved bulk sanitizer. Sanitizer top-up performed by the user.

LIFE+GUARD needs no installation. Cost effective with rapid deployment.

## BOX CONTENTS



- LIFE+GUARD Dispenser
- USB 2.0 Charging Cable
- Belt Clip
- Lanyard Clip with Lanyard
- User Manual
- Black and White colour options available



Belt Clip



Lanyard

**Marketed by:**

**Majac Medical Products Pty. Ltd.**

**Tel: 1300 138 578 | Int: +61 7 3265 6355**

**Email: [sales@majacmedical.com.au](mailto:sales@majacmedical.com.au)**

**[www.majacmedical.com.au](http://www.majacmedical.com.au)**

