SAPPHIRE [AMC 2130]

Suitable for Gynaecology and Obstetrics, with a variety of customisable options that allow you to get the functionality you need.

•••••

Electric adjustments of height and back rest.

Designed to reduce OH&S risks and minimise manual handling.



5 Year warranty on electrical components



warranty on frame and workmanship

Made



*Pictured paper roll holder and collection tray not included as standard. Available as optional customisations

Features

Movements powered by high quality electric actuators, reducing OH&S risks.

Supplied with foot control switch for hands free operation.

Supplied with either padded foot rests or under knee stirrups.

Low minimum height gives patients the ability to easily get on and off the couch, including wheelchair transfer.

Carefully considered base design allows easy side access for patients and practitioners.

Access under couch for patient lifting devices.

125mm lockable wheels allow effortless movement of the couch.

Basic specifications

Width: 610mm

Length: 1350mm

Back rest tilt: 75°

Height adjustment: 450mm -1000mm

Lifting capacity: 200kg

Wheels: 125mm lockable castors

Optional customisations

Removable foot section to convert couch to a longer, general purpose

examination table

Backward couch tilt (Trendelenburg) Stainless steel collection tray

with sliding arms and brackets, removable for cleaning

Extra width - 660mm

Battery pack

Hand switch control

Head cushion

Paper roll holder

Wide variety of fabric colours

Padded swivel arm rest set, suitable for Phlebotomy

Standard arm rest set

White LED light, attached to the couch or to a moveable stand

Central wheel locking



Majac Medical Products Pty Ltd


Ph: +617 3265 6355 Fax: +617 3865 2729 sales@majacmedical.com.au www.majacmedical.com.au




Also Available From:

SAPPHIRE (AMC 2130)

Top view



Side view



Optional customisations



Under knee stirrups



White LED light



Removable end section



Paper roll holder

Detailed specifications

Head section length: 870mm

Middle section length: 450mm

Foot section length: N/A

Vinyl specifications:

UV resistant

Flame resistant

Easy to clean

Sanitised to inhibit growth of bacteria

Complies with AS 1530.3 standards

Variety of colours available

Foam specifications: 50mm foam

Frame material and colour:

Powder coated white frame

Australian steel

Wheel system:

Individual 125mm brake castors made by Tente

Drive actuators:

Hi/lo: Dewert electric actuator system

Back rest: Dewert electric actuator system

Electrical specifications: 240V 50H 26VDC 105W

Electrical standards: IPX4 rating

Lift capacity: 200kg

Shipping weight: 70 - 85 kg

Shipping dimensions:

2130mm x 820mm x 700mm

Freight charges:

Not included in sale price.

Contact our sales representatives for a quote.

ARTG number: 208309





Majac Medical Products Pty Ltd

Ph: +617 3265 6355 Fax: +617 3865 2729 sales@majacmedical.com.au www.majacmedical.com.au



Also Available From: