



# SafeMask<sup>®</sup> Architect Pro<sup>™</sup>

Surgical N95 Respirator



MADE IN CANADA  
**BY THE MEDICAL  
MASK EXPERTS**

TGA - ARTG No. 380210


# Reliable respiratory protection that filters out over 95% of particulate aerosols free of oil

Adjustable nosepiece

Unique single-split comfort bands

Lightweight, breathable construction

High-quality materials and duck-bill shape



Made according to recognized biocompatibility, flammability and fluid resistance standards for use in healthcare facilities, the Medicom SafeMask® Architect Pro™ ensures reliable respiratory protection by filtering more than 95% of airborne particles free of oil. Its innovative “duck bill” shape provides optimal breathability with no exhalation valve required. Provides up to 8-hours of protection. Available in 3 sizes.

Medicom SafeMask Architect Pro N95 surgical respirator has full approval from the National Institute for Occupational Safety and Health (NIOSH) of the Centers for Disease Control and Prevention (CDC) under approval number TC-84A-8375. NIOSH certification means that the respirator has met the requirements of Title 42, Code of Federal Regulations, Part 84 (42 CFR 84).

TGA - ARTG No. 380210



#### Item codes:

**203214** Small      **203314** Medium      **203414** Large  
50 masks/box; 10 boxes/case



Majac Medical Products Pty Ltd

Ph: +617 3265 6355  
Fax: +617 3865 2729  
sales@majacmedical.com.au  
www.majacmedical.com.au



FS 540017

Also Available From:

