

Sontara

MediClean®

Low Lint Performance

Made from Sontara, new Mediclean towels blend the strength of polyester with the absorbency of cellulose. They provide a soft, highly absorbent towel with low lint and high wet strength properties, at a lower cost.

- Low lint performance
- Soft and highly absorbent soaking up to 7 times its own weight in fluids
- Superior dry and wet strength
- Steam stable
- Individual towels will not separate into layers
- Large enough to use double thickness
- Available in dispenser packs- 2 Sizes

Prefect for :

- Surface Cleaning
- Patient Care
- Instrument and Equipment drying



Code:

- 57-112 31.8x30cm (12 boxes/carton)
- 57-112/1 31.8x30cm (100 towels/ box)
- 57-113 35x60cm (12 boxes/carton)
- 57-113/1 35x60cm (100 towels/ box)