

Plasmapp Sterilization Product Catalog

Fast Low-Temperature Plasma Sterilizer “ 55°C ”

STERLINK®



Chamber : 14 Liter
3 Cycles Modes :
 - Fast : 7 min
 - Normal : 14 min
 - Full Load : 34 min
Class IIb, CE Certified
 * All cycles are Full Lumen Cycles

2019 New Model STERLINK® Mini



Chamber : 7 Liter
2 Cycle Mode :
 - Fast : 7 min
 - Full Load : 19 min
Class IIb, CE (June, 2019)
 * All cycles are Full Lumen Cycles

Specifications



Model	Dimensions				Weight [kg]	Chamber [liter]	Cassette (Sterilization cycle time)		
	Height	Width	Depth	Volume			STERPACK®	STERPACK® Plus	STERLOAD®
STERLINK®	437	433	614	116	68	14	○ (7 min)	○ (14 min)	○ (36 min)
STERLINK® Mini	330	274	440	40	19	7	○ (7 min)	X	○ (19 min)

Plasmapp Sterilization Product Catalog

▣ Sterilant Cassettes – H₂O₂ Included

STERPACK®



STERLINK® / STERLINK® Mini
7 min cycle
50pcs / Box

STERPACK® Plus



STERLINK®
14 min cycle
30pcs / Box

STERLOAD®



STERLINK® / STERLINK® Mini
36 min cycle / 19 min cycle
30pcs / Box


▣ STERPACK® Tray



In pouch aluminum tray for
STERPACK®



In pouch aluminum tray for
STERPACK® Plus



Chamber aluminum tray for
STERLOAD®

“ H₂O₂ Direct Injection Technology ”

“ World First Vacuum Sealed Pouch
Sterilization in 7 Minutes ”

“ Cycle Traceability System with
Individual Barcodes ”



Plasmapp Sterilization Product Catalog

▣ Sterilization Validation

Chemical Indicator



VH₂O₂ Class 1
CI Roll



VH₂O₂ Class 4
CI Strip

Biological Indicator



Compatibility :
3M, GKE, ASP

▣ Plasmapp VH₂O₂ Tyvek® Pouch Pre-Cut

TYVEK® 100



120pcs / Box
100mm X 400mm

TYVEK® 200



90pcs / Box
200mm X 400mm

TYVEK® 300



60pcs / Box
300mm X 400mm

▣ Accessories

CART



483x660x603mm

LABEL PRINTER



120x102x145mm

LABEL PAPER



Easy peel sticker

Temperature Control Auto-Rotary Sealer

STERSEAL® DE / STERSEAL® DE-V

One touch preset temperature button for STERPACK® mode.
Digital temperature controller for maximum adaptability



Validations

STERSEAL® DE

CE Directive 2006/42/EG
EMC Directive 2014/30/EU
WEEE Directive 2012/19/EU
RoHS 2011/65/EG

Harmonized Standards

EN ISO 12100/2010_11, 60204-1/2007_06
13857/2008_06, 61000-6-1/2007_10,
61000-6-3/2011_09

STERSEAL® DE-V

CE Directive 2006/42/EG
EMC Directive 2014/30/EU
WEEE Directive 2012/19/EU
RoHS 2011/65/EG

Harmonized Standards

EN ISO 12100/2010_11, 60204-1/2007_06
13857/2008_06, 61000-6-1/2007_10,
61000-6-3/2011_09

KRINKO/BfArM

EN ISO 11607-2:2006

ANSI/AAMI/ISO 11607-2:2006

DIN 58953-7:2010



Marketed by:
Majac Medical Products Pty. Ltd.
Tel: 1300 138 578 • Int: +61 7 3265 6355
Email: sales@majacmedical.com.au
www.majacmedical.com.au



FS 540017

Also available from:

