

# Apply & Wipe

Perfect Cleaning Partners

## ProDet®

Clinical Detergent

An environment cleaner formulated to solubilise and remove organic matter.

ProDet is "The Benchmark for Clinical Detergents". A powerful formula for all Healthcare requirements. Dynamic on Instruments, devices and surfaces.



- Multi-purpose Instruments & Surfaces
- Standards compliant
- Non-Hazardous
- Three year shelf life

### Product Features:

#### Excellent detergency with "rapid wetting"

Detergency is critical for removing bio-soil from the surface being cleaned. "Rapid wetting" is important for ultrasonic cleaning.

**Surface Cleaning is an essential** part of a complete hygiene protocol. Organic matter can inactivate many disinfectants and cleaning reduces the bio-burden.

**ProDet replaces a negative charge on all surfaces** this prevents re-deposition.

**Good Infection Control requires the removal of contaminates** such as bacteria and viruses with thorough cleaning then followed by disinfection, if necessary.

### Order information:

|                                 |   |
|---------------------------------|---|
| PD1- ProDet 1L                  | PDP1- 8ml pump for 1L container               |
| PD5- ProDet 5L                  | PDP5- 8ml pump for 5L container               |
| PD15- ProDet 15L (Tap included) | PDDS- 500ml dispenser Spray Top "non-misting" |
|                                 | PDDP- 500ml dispenser Pour Top                |

## Sontara MediClean®



### Low Lint Performance

Made from Sontara. The MediClean towels blend the strength of polyester with the absorbency of cellulose. They provide a soft, highly absorbent towel with low lint and high "wet strength" properties.

- Low lint performance.
- Soft and highly absorbent soaking up to 7 times its own weight in fluids.
- Superior dry and wet strengths. Steam stable.
- Individual towel will not separate into layers.
- Large enough to use double thickness.

### CODE:

57-112 31.8 x 33 cm  
(100 towels/box) (12 boxes/carton)  
57-113 35 x 60 cm  
(100 towels/box) (12 boxes/carton)

Apply ProDet to soften, lift & suspend bio-soils .....  
Wipe with MediClean to clean & dry

Link: ["Coronavirus \(COVID-19\) environmental cleaning and disinfection principles for health and residential care facilities"](#)

Marketed by:

Majac Medical Products Pty Ltd

Tel: 1300 138 578 \* Int: +617 3265 6355

Email: sales@majacmedical.com.au

[www.majacmedical.com.au](http://www.majacmedical.com.au)



FS 540017

Also available from:

