

# Apply & Wipe

"Healthcare facilities are concerned about cleaning sensitive materials with various chemical solutions".

## PlasClean®



**REMOVES THE BIO-BURDEN AND PREVENTS ENVIRONMENTAL STRESS CRACKING (ESC) ON PLASTIC OR VINYL SURFACES**

PlasClean® is an exceptionally environmentally friendly product as it is formulated from advanced surfactants and sequestering agents that are made using renewable biological resources.

### Product Feature:

#### Clean & Protect:

- Dental Chairs,
- Podiatry Chairs,
- Treatment/Exam Beds,
- Operation Tables,
- Examination Lights,
- Medical Devices.

#### PlasClean® is also ideal for:

- Connectors, lines, and tubes,
- Ophthalmic equipment,
- Training manikins,
- Respiratory apparatus, such as resuscitators, masks, pumps, and nebulizers.

#### Order information:

PLC5- PlasClean 5L

PLC500T- 500ML dispenser (spray top). Empty.

PLC500P- 500ML dispenser (pour top). Empty.

## Sontara MediClean®



### Low Lint Performance

Made from Sontara. The MediClean towels blend the strength of polyester with the absorbency of cellulose. They provide a soft, highly absorbent towel with low lint and high "wet strength" properties.

- Low lint performance.
- Soft and highly absorbent soaking up to 7 times its own weight in fluids.
- Superior dry and wet strengths. Steam stable.
- Individual towels will not separate into layers.
- Large enough to use double thickness.

#### CODE:

57-112 31.8 x 33 cm

(100 towels/box) (12 boxes/carton)

57-113 35 x 60 cm

(100 towels/box) (12 boxes/carton)

**Apply PlasClean to soften, lift & suspend bio-soils .....  
Wipe with MediClean to clean & dry**

Link "[All about Environmental Stress Cracking \(ESC\). PlasClean](#)"

Marketed by:

Majac Medical Products Pty Ltd

Tel: 1300 138 578 \* Int: +617 3265 6355

Email: sales@majacmedical.com.au

[www.majacmedical.com.au](http://www.majacmedical.com.au)



Also available from: