

## Bowie Dick Test Pack



SAP / Monitoring

Once a day ..... since 1963

REF: 102.001.0020 Unit: 20 pieces

Used for daily monitoring of pre-vacuum steam sterilizer. Based upon the Bowie Dick test result the steam sterilizer is released for usage during the day.

Characteristics:

- ◆ Fulfill requirements of EN 285 and ISO 17665
- ◆ Conform to ISO 11140 - 1 class 2
- ◆ Independently tested to show performance similar to original pack
- ◆ Clear yellow - black color change
- ◆ Easy to interpret
- ◆ Easy to document
- ◆ Non-Toxic / Lead Free

The Bowie Dick test detect failures like air leaks, insufficient vacuum, inadequate steam penetration.

