


AMC 2120 Topaz

- Suitable for a range of disciplines and customisable to your needs
- Electric adjustments of height, backrest, independent leg sections and backward tilt (option)
- Designed to reduce OH&S risks and minimise manual handling
- Complies with RACGP guidelines on examination couches

Specifications

Safe working load: 200kg

Colours: ●●●●●●●●●●●●●●●●



AUSTRALIAN OWNED & MADE



10 YEAR WARRANTY
On frame and workmanship



WE CUSTOMISE
Built to suit your needs

Configure

- Select the accessories that suit your practice

Independent leg sections

- Electrically powered footrest height adjustment

Safety

- Mechanically and load cycle tested
- 125mm individual locking castors
- Electric locking castors (option)
- Automatic battery backup during power outages (option)

High-quality vinyl

- Multiple colours
- Flame- and UV-resistant
- Antimicrobial and antibacterial
- Highest level of stain resistance
- Easy to clean

Fully adjustable


- Adjust the position to work safely and comfortably
- High resolution height and backrest adjustment
- Performance electric actuators

Total control

- Easy to use controls
- Ergonomic design
- Handset and/or footswitch options

Maximum comfort

- Compact frame design to work closer to the patient
- Designed for patient and practitioner comfort
- Low minimal height for easy patient access
- Under couch access for patient lifting devices



Also available from:



Marketed by:
Majac Medical Products Pty. Ltd.
Tel: 1300 138 578 • Int: +61 7 3265 6355
Email: sales@majacmedical.com.au
www.majacmedical.com.au



FS 540017