



**Pacific Medical**  
Australia

## 3 Section Medical Treatment Couch with Side Rails

### Features:


- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

### Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m<sup>3</sup>
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

### ORDER CODES:

SRET33NB – Navy Blue



### 3 Section Medical Treatment Couch with Side Rails Demo Video

**MAINTENANCE:** Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.



Marketed by:  
**Majac Medical Products Pty. Ltd.**  
Tel: 1300 138 578 • Int: +61 7 3265 6355 Email:  
sales@majacmedical.com.au  
[www.majacmedical.com.au](http://www.majacmedical.com.au)



FS 540017

Also available from: