



Pacific Medical
Australia

Beauty Day Spa Couch

Features:

- Electric Hi Lo operation by foot switch
- Luxurious soft and comfortable head cushion
- Contoured comfortable table top
- Sleek design with rounded frame legs for safety
- Removeable padded wing armrests
- Padded prone arm rest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 186cm
- Width: 62cm
- Backrest length: 58cm
- Mid section length: 45cm
- Legrest length: 41cm
- Long section length: 121cm
- Adjustable backrest through gas operated strut to +85 Degrees
- Adjustable footrest through gas operated strut from +5 to -75 Degrees
- Height range from 57cm to 98cm in approx 30 seconds
- Padding thickness 7.5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

- BTWH - White
- BTNB - Navy Blue

[Beauty Day Spa Couch Demo Video](#)

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.



Marketed by:
Majac Medical Products Pty. Ltd.
Tel: 1300 138 578 • Int: +61 7 3265 6355 Email:
sales@majacmedical.com.au
www.majacmedical.com.au



FS 540017

Also available from: